
F. BOSCH
MEDIZINTECHNIK

Model: **Provac**

PORTABLE PHLEGM SUCTION MACHINE
ELEKTRİKLİ ASPIRATÖR CİHAZI

1 Product Usage	3
2 Structural characteristics and working principle	3
3 Main Technical Specifications	4
4 Standard operating conditions	4
5. Commissioning	4
5.1 Opening Check	4
5.2 Tubing Connections	4
5.3 Power Connections	6
5.4 Tubing calibration	6
5.5 Pressure regulation	6
5.6 Inspection and testing of overflow component	7
5.7 Shutting off	7
5.8 Cleaning and Hygiene.....	8
6 Safety symbols and their meanings used on this device	8
7 Operations and Maintenance	9
7.1 Operation and Maintenance	9
7.2 Change of air filter	9
7.3 change of fuse	10
7.4 Maintenance after use	10
7.5 Periodic Safety Checks	11
8 Common Troubleshooting	11
9 Additional Notes	12
9.1 Restriction on transportation and storage conditions	12
9.2 Electrical Theory Flow Chart	12
9.3 After sales service	12
9.4 Accessories list	13
10. Guidance and manufacture's declaration-electromagnetic emissions-for all EQUIPMENT and SYSTEMS	13
11.Manufacturer.....	14

Portable Suction Device Operation Manual

1. Product Usage

The Phlegm Suction Unit is a portable, AC-powered device intended to be used to remove Phlegm from a patient's airway. Unfit thoracic drainage and gynecological drainage.


The device may be used during surgery in the operating room or at the patient's bedside.

2. Structural characteristics and working principle

The main structure of suction device includes vacuum pumps, vacuum gauge, pressure regulating valve, air filter, and the reservoir bottles.

It used oil-free lubrication pump to stop environmental pollution by oil mist. It has low noise. Reservoir bottle is hidden in a submerged novel structure, and it has an overall plastic enclosure design. The operation of the device will not produce positive pressure to ensure reliability and safety. The vacuum conditioning system can be adjusted by non-polar voltage as needed. It is small size, light weight and easy to carry, and is suitable for a variety of emergency situations and out-clinics in particular.

Operation Theory flow chart

3. Main Technical Specifications

1. Power: AC 230V, 50Hz	2. Input current: 0.7A
3. Max negative pressure value : 87 Kpa (0.87 bar)	4. pressure adjustment range: 0.01 MPa ~ Max Pressure Value Limit
5. Pumping rate: ≤ 20 L/MIN	6. Fuse: AC250V, F2.0 AL, $\Phi 5 \times 20$
7. Reservoir bottle: 1000mL, One Piece	8. Noise: ≤ 60 dB(A)
9. Vacuum gauge: 0~100KPa, 2.5 grade	10. Equipment Type: Medium vacuum / low flow
11. Working mode: Short-time working	12. Electrical requirements: is a class II ,type BF
Note: It is not suitable for use in environment with flammable and explosive gases; Electromagnetic interference source: the minimum separation distance of ≥ 1.1 m when it is close to power lines and electrical equipment with the voltage less than 500V	

4. Standard operating conditions

Temperature: $+5^{\circ}\text{C} \sim +40^{\circ}\text{C}$	Relative Humidity Range: $\leq 93\%$
Atmospheric pressure range: 500hPa~1060hPa	Transport and storage temperature range: $-25^{\circ}\text{C} \sim +70^{\circ}\text{C}$
Note: When the device is stored and transported in temperature below 5°C , the device should be placed in normal working temperature environment for more than 4 hours before use.	


5. Commissioning

5.1 Opening Check

Before the device is installed and commissioned, the user should check whether the product Appearance look good, and whether type and quantity of accessories are consistent to the accessories list attached in the Manual. Please contact the supplier or manufacturer promptly if there is any shortage or defect.

5.2 Tubing Connections

(Please refer to tubing connection diagram, and note suction tube is not connected)
Note: Before installation, apply small amount of distilled water on the mouth
section where cap is pressured into the bottle to improve compaction and increase
the seal tightness.

- 1) Reservoir bottle 2) suction soft tube 3) air filter 4) blue mark
5) suction catheter

Tubing Connection Diagram

Tubing Connection Diagram

Accessories(1)

Tubing Connection Diagram(1)

5.3 Power Connections

Connect the power plug to power supply, and the power supply is connected when power indicator is on.

Note: The power cut-off device of this product is the power plug.

5.4 Tubing calibration

Tightened pressure regulating valve clockwise, block aspiratory port with fingers or rubber tip of dropper, or fold and pinch the suction catheter;

When the suction switch is on and machine is operated without abnormal sound, and the vacuum gauge pointer rise sharply to pressure value limit, and meter pointer return to below 10KPa when suction moth is released, the tubing is connected correctly.

Note 1: Suction catheter is not provided by us , Please use suction catheter with CE mark through your provider.

Note 2: If the suction tube is blocked, please clear the tube by following method: Fold the suction catheter into a V-type (reservoir bottle should be empty) for the pressure to reach its maximum value, then release the catheter quickly, repeat the same action several times to force the suction tube to flow.

5.5 Pressure regulation

Block suction mouth, turn on the suction device, and adjust the pressure regulating valve, the reading on vacuum gauge should fluctuate between 0.02MPa~Pressure value Limit.

During clinical use, the pressure regulating valve is used to control the suction pressure required.

Turn pressure regulating valve clockwise to increase pressure. The pressure must drop below 0.02MPa before turn off.

5.6 Inspection and testing of overflow component

Open the bottle, clean the valve port, and press on the rubber valves on flat floater. The valve should not have warping, tearing and other defects, and well connected with the floater. Floater in the float planes should move freely and without unimpeded rolling condition.

Hold up bottle cap to allow the floater a vertical contact with the liquid surface, the drop cap slowly, the floater should be able to float in the float planes.

Tighten the cap, connect the soft suction catheter to suction mouth, tighten control valve to operate the suction device. Extend suction catheter into a bucket of clear water or in a simulated normal usage environment, to suck the liquid into the reservoir bottle which has an overflow component, rise of liquid level will lead floater to float until the floater closes off at the valve port and stop the device automatically. The

final liquid level will change subject to different suction methods.

Loosen the valve, and turn off the suction switch, open the bottle, and empty reservoir bottle, when the cap is tightened again, the floater should be at the bottom of the float planes, the valve port was open. When the above readings are observed, the overflow component is working order and can be used in operation.

Note 1: If when the overflow component is shut off and the liquid level continues to rise, there might be two reasons:

- (1) Due to the pressure remained in the reservoir bottle;
- (2) Valve port is not completely closed.

In case (1) pull off the suction catheter away from liquid area of suction and put the catheter back into the liquid area of suction, the liquid level in reservoir bottle should stop rising; in case (2) the liquid level will continue to rise in which case, the situation should be carefully observed, when the reservoir bottle is near full, the suction catheter should be pull off from the liquid area of suction immediately, turn off the suction device, stop suction and find out the reasons for valve port failure.

Note 2: Valve port is closed off by the floater, and suction stopped. However, due to the pressure inside the tube, the floater may be stuck on the valve port. In this case, one should spin loose control valve or turn off suction devices to release pressure inside the tubes, and the floater should fall from the valve port by gravity. (Pulling the floater off by hand is prohibited to prevent the rubber valve of the floater being torn off. If there is any mucus on the floater, the floater should be thoroughly cleaned before it can be used again.)

Note 3: The bottle may be uncapped only after the device is shutdown and pressure is released.

Note 4: The suction device must not be used when overflow component is removed.

5.7 Shutting off

After installation, commissioning or usage, the suction device must be turned off before unplug the power plug to cut off power.

5.8 Cleaning and Hygiene:

5.8.1 Basic information:

The Phlegm Suction Machine is for patient care . It may share use with others, but the Suction catheter single use only. Except the Phlegm Suction Machine and air filter, the accessories can be washed and sanitized.

.

The accessories life is 3 years , the air filter need to change when turn yellow..

5.8.2 First use and once daily:

- . Disassemble all accessories to pieces and flush each piece by tap water. (All accessories mean the all parts except the Phlegm Suction Machine and air filter).
- . Connected the tube to the machine;.

5.8.3 After each use:

- . Clean the accessories. Wash by warm soapy water and rinse by clean water.

5.8.4 Keeping care:

- . To care the accessories a good hygiene, suggest to put the accessories in a container or a plastic bag with seal until

5.8.5 About Phlegm Suction Machine:

- . Do not pour water onto Phlegm Suction Machine. Do not immerse the Phlegm Suction Machine into water. If there are some water drops onto Phlegm Suction Machine, use a wet towel to clean surfaces.

6. Safety symbols and their meanings used on this device

Symbol	Meaning	Symbol	Meaning
~	AC		Consult Accompanying documents
	Class II device		Type BF
O	OFF	I	ON

7 .Operations and Maintenance

7.1 Operation and Maintenance

Prior to operation, check the suction device according to the installation and commissioning procedures to ensure it is in good working order, and once sterilized suction catheter and tube is connected, the device may be put into use.

Note: please read the packaging description before using the attached suction tube.

During operation, pressure regulating valve is used to adjust the pressure to desired level, and turned the switch on or off as per requirement. Always pay attention to liquid level in the reservoir bottle. When the liquid level rises to the calibration capacity of the reservoir bottle (the device remain operational when tilted 10 degrees), suction operation should stop to empty and clean reservoir bottle before use again, otherwise the liquid level will drive the floater up to the level that closes off valve port, forcing suction operation to stop automatically.

Note: When liquid level continues to rise after overflow component is shut off, Please refer to inspection and testing of overflow component procedure for handling.

Emergency measures in use:

- 1) In case when purulent sputum or mucus obstructs suction tube, the pressure regulating valve should be loosen promptly to release pressure, and suction tube should be replaced before used again;
- 2) In case where the suction tube difficult to pull out after suction or sucked on human tissues, please apply the method above by loosening the pressure regulating valve.

Operating time must not exceed 30 minute

Suction volume must not exceed the maximum liquid level warning mark

Tightened cap and tube joints when pressure is insufficient

Note 1: Prior to suction, bend the suction catheter into a V shape to boost that pressure to the required range on turning on the device, then insert the suction tube into patients Phlegm site, and then release the suction catheter to suck for better result.

Note 2: The specification and size of suction tube is determined by the medical staff according to clinical requirements

Note 3: During operation, if there is presence of bubble in the reservoir bottle, liquid level of the reservoir battle should be monitored to prevent overflow device failure.

Note 4: The suction devices should be used under the guidance of medical personnel, and in strict accordance to operation procedures outline in the Manual, and in case of doubt please contact the supplier or manufacturer.

7.2 Change of air filter

When air filters sucked or filled with dust would result in colour change of membrane from light to dark, and also caused significant reduction and even disappearance of inhalation capacity at the tube mouth, and pressure on the vacuum gauge will continue to rise to more than 0.04MPa. In such case, please replace the air filter replacement from our company.

Note 1: When the overflow component is shut off or the suction tube is blocked, suction force may also decrease or disappear, and pressure values would also increased. Please refer to “Common Troubleshooting”.

Note 2: The air filter should be changed regularly and destroyed in a centralized method.

7.3 Change of fuse

Fuse is located in the bottom of base, when replacement is required, cut off the power supply first, use flat or cross screwdriver counter-clockwise to open the base to replace.

7.4 Maintenance after use

Prior to shut down, Wash in a hot water/dishwashing detergent solution and rinse with clean, hot tap water. Then wash with a commercial (bacterial-germicidal) disinfectant, follow disinfectant manufacture's recommended instructions and dilution rates carefully.

It is recommended to suck small amount of clean water using suction tube to clean the tube internal wall. After shut down, empty reservoir bottle, use soft brush or cloth to remove dirt on bottle and the cap, then wash again with water (including the overflow component, seals and various tubes). If necessary, unscrew the overflow component, separate float planes and floater for a through cleaning (Note: rubber valve shall not be separated and floater).

Used suction tube should be cleaned by saline to remove remnants of purulent sputum and mucus within the tube. If the suction tube is blocked, it should be replaced.

(Suction catheter is not provided by us , Please use suction catheter with CE mark through your provider.)

Reservoir bottles, caps and various tubes should be soaked and disinfect in three parts hot water to one part white vinegar for 30 minutes or use a medical bacterial-germicidal disinfectant available through your provider.

Note: Reservoir bottle is made of plastic, avoid sharp items when clean or in use, and also prevent it from falling

The outer surface should be wiped by wet cloth pre-soaked with disinfectant. To avoid the infiltration of liquid, do not wipe on crack, words and artworks
Devices not in use should be placed in dry and clean place, and operate the device on a regular basis (typically every six months).

Note: Before the next use of the device, overflow component and other tubing should be fitted according to tube connection method.

7.5 Periodic Safety Checks

The following safety checks should be performed at least every 24 months by a qualified person who has adequate training, knowledge, and practical experience to perform these tests.

- Inspect the equipment and accessories for mechanical and functional damage.
- Inspect the safety relevant labels for legibility.
- Inspect the fuse to verify compliance with rated current and breaking characteristics.
- Verify that the device functions properly as described in the instructions for use.
- The leakage current should never exceed the limit defined in EN/IEC 60601-1. The data should be recorded in an equipment log. If the device is not functioning properly or fails any of the above tests, the device has to be repaired.

8 .Common Troubleshooting

Ref	Failure	Cause analysis	Remedy	Note
1	Negative Pressure limit is less than 0.07MPa	1)mouth leak 2)tube junction leak 3)regulating valve loose or open	1)clean bottle mouth, tighten or replace the cork, seals or joints 2)re-plugged the junction 3)tighten control valve	1)maintenance of cabin parts should be conducted by qualified personnel; 2)split tube should be replaced
2	Negative value is greater than 0.04MPa, but the suction power at mouth of the suction tube is significantly reduced or disappeared	1)Overflow component is shut off 2)Tubing is blocked 3)Air filter is blocked	1)turn off the device, turn the regulating valve counter-clockwise to release the pressure in the tube, then tighten the valve 2)Dredge, clean or replace the tube 3)Replace the air filter from Our company	1)empty the reservoir bottle promptly 2)blue mark on the air filter is the inhaling mouth

3	Power supply voltage is normal, but indicator light is not on	1)power plug is loosen 2)fuse is fused 3)indicator light is damaged	1)repair or replace power plug 2)replace fuse 3)replace indicator light	1)specification:AC250V, F1.5AL, Φ5X20
4	Fuse is fused	1)voltage oversupplied 2)internal line fault 3)rolling resistance of pump and resulted in current increase	1)adjust voltage 2)check and repair line fault 3)check pump and motor	Must be conducted by qualified technician (refer to electrical Theory flow Chart)


Note: If pump fails, demolition and repair should be operated by qualified personnel, please contact the manufacturer if necessary (power should be cut off before inspection of lines or opening up the box)

9 .Additional Notes

9.1 Restriction on transportation and storage conditions

Portable suction device should be stored well-ventilated room where relative humidity is no more than 80% and without corrosive gases, and should avoid sever vibration or shock.

9.2 Electrical Theory Flow Chart

Electrical maintenance should be conducted by qualified personnel

9.3 After sales service

The company will accept return, replace or repair for product with quality issues that was not caused by artificial factor and occurred within one week from the date of sales. The company will provide maintenance free of charge for product with quality issues which occurred within one and a half from the date of manufacturing (one year from use commence date) under standard use and storage conditions. For products with quality issues which occurred after one and a half from the date of manufacturing (one year from use commence date), the customer should take invoice and warranty

card to our after-sales service department, representative office or dealer office, and the company will provide parts for maintenance at a reasonable charge. Where user is not able to provide invoice, the company will extend the free warrantee period by one month calculated based on device reference or date of manufacturing.

Warranty does not extend to the following circumstances: 1. machine damage or deformation caused by collision; 2. machine exposed to water or rain; 3. infiltration of water, blood or sputum, or suction of viscous liquid into the pump by customer and resulted in machine failure; 4. vulnerable and consumables items including: air filters, suction tube and fuse.

When product and accessories are reaching the end of their expiry dates, these products and accessories should be disinfected then be delivered to qualified companies or individuals for recycling.

9.4 Accessories list

Suction tubing(length1.3m)	1 PCS	Intermediate tubing (length 0.2m)	2 PCS
Manual	1 PCS	Air filter	3 PCS
Reservoir bottle	1PCS	suction catheter(12F/8F)	2PCS
Catheter connector	1PCS		

10.Guidance and manufacture's declaration-electromagnetic emissions-for all EQUIPMENT and SYSTEMS

Guidance and manufacture's declaration-electromagnetic emission		
The BR-SM181 is intended for use in the electromagnetic environment specified below.The customer of the user of the BR-SM181 should assure that it is used in such and environment.		
Emission test	Compliance	Electromagnetic environment-guidance
RF emission CISPR 11	Group 1	The BR-SM181 uses RF energy only for its internal function.Therefore,its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class B	
Harmonic emissions IEC61000-3-2	Complies	
Voltage fluctuatons/flicker emission IEC61000-3-3	Complies	The BR-SM181 is suitable for us in allestablishments,including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes

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